

# cuccio PRO™



Cures 5 fingers.  
LCD control panel with On/Off button and timer. Removable tray.  
Power output 18 W.  
Wavelength 395nm to 405nm



Cures 5 fingers. rechargeable  
Digital display  
Comfort Cure mode. Removable tray.  
Power output 18 W.  
Wavelength 395nm to 405nm








Cures 5 fingers.  
LCD control panel with On/Off button and timer. Removable tray.  
Power output 18 W.  
Wavelength 395nm to 405nm



Cures 4 fingers.  
On/Off button. Cures Most gel polish, Top Coat, Clear hard gels.  
Power output 6W.  
Wavelength 395nm to 405nm



Cures 4 fingers.  
On/Off button. Cures Most gel polish, Top Coat, Clear hard gels.  
Power output 6W. USB cable  
Wavelength 380nm to 410nm



# cuccio PRO™



Cures 5 fingers.  
LCD control panel with On/Off button and timer. Removable tray.  
Power output 18 W.  
Wavelength 395nm to 405nm



Cures 5 fingers. rechargeable  
Digital display  
Comfort Cure mode. Removable tray.  
Power output 18 W.  
Wavelength 395nm to 405nm



Cures 5 fingers.  
LCD control panel with On/Off button and timer. Removable tray.  
Power output 18 W.  
Wavelength 395nm to 405nm



Cures 4 fingers.  
On/Off button. Cures Most gel polish, Top Coat, Clear hard gels.  
Power output 6W.  
Wavelength 395nm to 405nm



Cures 4 fingers.  
On/Off button. Cures Most gel polish, Top Coat, Clear hard gels.  
Power output 6W. USB cable  
Wavelength 380nm to 410nm

